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
---


## *HIGH SPEED BEARINGS*

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*A dental bearing must be able to withstand the most extreme conditions and is, at the same time, one of the most important components for the performance of a dental instrument. Tecnomed chooses only the best bearings in the dental industry market, with both ceramic and steel balls. These are precise and reliable, thanks to the use of the highest quality materials that guarantee a long product life. Minimum noise emission is ensured even at high speeds. High-end bearings are sold in a sealed package and marked with their country of origin.*

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*For decades, Myonic has been supplying customers around the world with its innovative and reliable dental bearings, gaining an exceptional position in the market. Rotational speeds of up to 500,000 rpm, vibrations, overheated steam or various cleaning tools, as well as the penetration of dirt, such as blood and tooth powder residue, are the biggest challenges for these special bearings.*

*In addition to a unique design, these challenges are met with a systematic selection of materials and lubricants. The result is a product tailored to this application that meets all the desired requirements over time.*

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*Searching for high-quality products on the market at competitive prices is our job. For over 10 years, Tecnomed has also been offering bearings under its own brand, thanks to a partnership with an important North American company. Our dental bearings boast of a high quality and precision. They are ideal for use in dental handpieces, medical devices, and surgical instruments. With a grade 3 ball and an accuracy that is higher than an ABEC 7 bearing, they are very reliable. The bearings are sold with an industrial packaging.*

myonic

optimyn

TecnoMed  
ITALIA Group

## FENOLICI O TORLON? FENOLIC OR TORLON?

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*The ball bearing cage is designed to keep the ball in the ball bearing separate from one another, around the partial rim of the bearing. To be able to find the optimal solution for any ball bearing, Myonic has developed several models of ball bearing cages. There is no single cage that meets all the requirements, therefore, the choice of an optimal cage must be made according to specific needs. There are two types of cages that perform better in the dental sector: phenol resin cages and torlon cages.*

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
*Phenol resin cages, reinforced with fabric (brown color), are very light and can, at the same time, resist acceleration forces and very high centrifugal forces.*

*PRO: They retain lubrication fluid longer*

*CON: They suffer most during the sterilization cycle*

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*The use of torlon (black color) for cages is increasingly widespread in applications involving high speeds, chemical corrosion, or high temperatures.*

*PRO: The exceptional properties of torlon are due to the optimal combination of strength, flexibility, compatibility with high operating temperatures, good workability, and high chemical and wear resistance.*

*CON: Since it is not a porous material, it does not retain lubricating fluid for long, therefore, the handpiece must be properly lubricated.*

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*Ceramic ball bearings are aesthetically similar to those with steel balls. The only substantial difference is that the balls inside the bearing are ceramic rather than steel.*

*Plates and cups are made from ceramic. So, what sense does using this material have? Actually, the same type of ceramic is not used for ball bearings. So-called ceramic materials, and, more precisely, Silicon Nitride (Si<sub>3</sub>N<sub>4</sub>), are used.*

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- *These bearings are quieter: ceramic balls do not become deformed because they are harder and, even after many hours of work, the bearings do not wear out and continue with their original level of silence. Instead, steel balls become deformed over time and, when worn out, both friction and noise increase.*
- *They also respond better to high temperatures. Therefore, during a sterilization cycle, they do not become deformed, they maintain their shape, and the noise emitted does not increase. Steel balls, on the other hand, dilate and undergo deformation at high temperatures. Consequently, both friction and noise increase.*
- *They are lighter than steel ones, even if slightly. This helps lighten the tool, making it more manageable.*
- *They do not oxidize because ceramic materials are not sensitive to water, steam, and humidity. On the contrary, steel balls are subjected to oxidation.*
- *They better dissipate the heat that is generated during the working phase, ensuring a longer life and a low handpiece temperature.*

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
*The difference between radial and angular contact ball bearings lies in the type of internal sealing ring used, which defines the axial and radial load bearing support.*

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
*Radial contact ball bearings are those where the load force to be supported is perpendicular to the rotation axis. These bearings are particularly versatile. They are also suitable for use at high and very high speeds, they can withstand axial and radial loads in both directions, and require little maintenance. Easy to handle, radial contact ball bearings are preferred by technical service centers, which is why Tecnomed Italia has chosen to offer only radial contact ball bearings.*

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*Angular (oblique) contact ball bearings have inner and outer raceways displaced relative to each other in the direction of the bearing axis. Given this conformation, these bearings can bear combined loads, i.e. acting simultaneously in radial and axial directions. Technically, angular contact ball bearings have a better technical performance but are sensitive to errors during assembly on the rotor, since they can easily come apart.*

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*The ABEC scale is a tolerance classification system used for precision bearings. ABEC is the acronym of the American organization that gives its name to this system: Anular Bearing Engineering Committee. This organization is part of the American Bearing Manufacturers Association.*

*Bearings classified with the ABEC system are called "precision bearings" and are marked with an odd number from 1 to 11. The highest value corresponds to the best precision standards:*

**Higher number = lower tolerance = more expensive bearings**

*All bearings for dental turbines on the market are ABEC 7 class bearings. ABEC class 9 dental bearings are rarely found. The latter are much more expensive.*

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Buon lavoro!

*What Tecnomed Italia offers are bearings that our technical service uses for daily repairs: Myonic bearings in a phenol cage with ceramic or steel balls and Tecnomed Italia ball bearings with a torlon cage and steel and ceramic balls. Next to each code, we indicate the brand, type of cage, and type of balls. Now that you are familiar with everything, you can select the most suitable bearing for your needs!*

*Have a nice day!*



**Tecnomed Italia Staff**  
**Technical service department**

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TECNOMED ITALIA Order code	Brand	Cage	Balls
CN02A	Tecnomed Italia	torlon	steel
CN22BC	Tecnomed Italia	torlon	ceramic
CN22MC	Myonic	fenolic	steel
CB22BC	Myonic	fenolic	ceramic

CROSS REFERENCE BRAND CODE				
Brand	Tecnomed Italia	Myonic	Timken	Barden
Steel bearing code	CN02A	CN22MC	DR02B2L	SR144STAY134K3C10GJ2
Ceramic bearing code	CN22BC	CB22BC	DR02A2L-801	CSR144STAY134K3C10GJ

CN02A / CN22BC / CN22MC / CB22BC | 6,35 x 3,175 x 2,380 mm smooth  
Applications


HIGH SPEED BEARINGS

Handpiece brand	Handpiece model	Front Bearing (burs side)	Rear Bearing (push side)
CHAMPION		•	
CHIRADENT	3200	•	•
	3400	•	•
FARO	S380	•	•
KMD	COLIBRI	•	•
KAVO	8000/B (rotor compatible TMI)		•
	7000/B (rotor compatible TMI)	•	•
IVORY-DENT	T501/L-SM502/4M		•
LARES	330E	•	•
	557 SMALL HEAD	•	
	757 LARGE HEAD	•	
MIDWEST	QUIET AIR LEVER/PUSH BUTTON		•
	QUIET AIR STANDARD		•
	TRADITION L/PB CONVERSION		•
	TRADITION STANDARD		•
MTC	105	•	•
NSK	CH-QD		•
	MACH LITE ML		•
	MACH MC		•
	MC-QD		•
	N75		•
	NL-75		•
	P&S	•	•
	PANA-AIR MINI PUSH BUTTON	•	•
	PANA-AIR MINI SCREW TYPE	•	•
	PANA-AIR STANDARD PUSH BUTTON	•	•
	PANA-AIR STANDARD SCREW TYPE	•	•

NSK	PANA-AIR TORQUE PUSH BUTTON	•	•
	PANA-AIR TORQUE SCREW TYPE	•	•
	PTL-II		•
	SUPER GRADE	•	•
	SUPER GRADE PB	•	•
	TRIPLE SPRAY/KINETIC QUANTUM		•
	VIP	•	•
TECNOMED ITALIA	HCP115/116/117/118	•	•
TKD	MINI MASTER	•	•
W & H (ADEC)	200 MINI FRICTION GRIP	•	•
	200 MINI PUSH BUTTON	•	•
	200 STANDARD FRICTION GRIP	•	•
	200 STANDARD PUSH BUTTON	•	•
	300 MINI FRICTION GRIP	•	•
	300 MINI PUSH BUTTON	•	•
	300 STANDARD FRICTION GRIP	•	•
	300 STANDARD PUSH BUTTON	•	•
	395-398-695-795-895	•	•
	TK94 L/LM	•	•
SIEMENS	3000	•	•
SIRONA	T2/T3 BOOST (Push with 3 cuts and rotor compatible TMI)	•	•
	T2 CONTROL (Push with 3 cuts and rotor compatible TMI)	•	•
	T2 CONTROL ENJOY	•	•
YOSHIDA	H-Q (FG)	•	•
	H-QP (PB) (HIGH TORQUE)	•	•
	SL-Q (FG)	•	•
	SL-QP (PB) (MINI HEAD)	•	•

CU8B / CT88C / CUT8T / CB89C | 7,94 x 3,175 x 3,571 mm smooth

HIGH SPEED BEARINGS

Size (mm)	
D	7.94
d	3.175
B	3.571

TECNOMED ITALIA Order code	Brand	Cage	Balls
CU8B	Tecnomed Italia	torlon	steel
CT88C	Tecnomed Italia	torlon	ceramic
CUT8T	Myonic	fenolic	steel
CB89C	Myonic	fenolic	ceramic

CROSS REFERENCE BRAND CODE				
Brand	Tecnomed Italia	Myonic	Timken	Barden
Steel bearing code	CU8B	CUT8T	DR09B2G	SR2-5STAY36C10GJ295E
Ceramic bearing code	CT88C	CB89C	DR09B2G-801	CSR2-5STAY36C10 0-10


CU8B / CT88C / CUT8T / CB89C | 7,94 x 3,175 x 3,571 mm smooth  
Applications

HIGH SPEED BEARINGS

Handpiece brand	Handpiece model	Front Bearing (burs side)	Rear Bearing (push side)
BIEN AIR	BLACK PEARL ECO	•	•
	BORA/L - S36L	•	•
	EOLIA	•	•
	BORALINA	•	•
	CROMA (STANDARD HEAD)	•	•
	LAB	•	•
BUFFALO	220 LAB	•	•
KMD	MAXTORQUE	•	•
NEYTECH	GRAND HURRICANE	•	•
	HURRICANE	•	•
	QC-700	•	•
STAR	CONCENTRIC	•	•
	FUTURA	•	•
TECNOMED ITALIA	HCP01P	•	•
TKD	MASTER	•	•

CU61S / CU61CC | 6,35 x 3,175 x 2,38 mm groove

HIGH SPEED BEARINGS

Size (mm)	
D	6.350
D'	6.00
d	3.175
B	2.38
B'	1.10

TECNOMED ITALIA Order code	Brand	Cage	Balls
CU61S	Tecnomed Italia	torlon	steel
CU61CC	Tecnomed Italia	torlon	ceramic

CROSS REFERENCE BRAND CODE				
Brand	Tecnomed Italia	Myonic	Timken	Barden
Steel bearing code	CU61S	/	DR21B2L	/
Ceramic bearing code	CU61CC	/	DR21B2L-801	/


CU61S / CU61CC | 6,35 x 3,175 x 2,38 mm groove  
Applications

HIGH SPEED BEARINGS

Handpiece brand	Handpiece model	Front Bearing (burs side)	Rear Bearing (push side)
DABI ATLANTE	MS350	•	•
SIEMENS	4000 M	•	•
	TM1	•	•

CU757L | 6,35 x 3,175 x 1,91 mm smooth

HIGH SPEED BEARINGS

Size (mm)	
D	6.35
d	3.175
B	1.91

TECNOMED ITALIA Order code	Brand	Cage	Balls
CU757L	Tecnomed Italia	torlon	steel

CROSS REFERENCE BRAND CODE				
Brand	Tecnomed Italia	Myonic	Timken	Barden
Steel bearing code	CU757L	/	DR74A1L	/

CU757L | 6,35 x 3,175 x 1,91 mm smooth  
Applications


HIGH SPEED BEARINGS

Handpiece brand	Handpiece model	Front Bearing (burs side)	Rear Bearing (push side)
LARES	557 SMALL HEAD	•	•


CU7C / CU77C | 6,35 x 3,175 x 2,78 mm flaged

HIGH SPEED BEARINGS

Size (mm)	
D	6.35
d	3.175
D'	7.5
B	2.78
C	0.80

TECNOMED ITALIA Order code	Brand	Cage	Balls
CU7C	Tecnomed Italia	torlon	steel
CU77C	Tecnomed Italia	torlon	ceramic

CROSS REFERENCE BRAND CODE				
Brand	Tecnomed Italia	Myonic	Timken	Barden
Steel bearing code	CU7C	/	DR01B2L	SFR144RSTAY167K4C10G
Ceramic bearing code	CU77CC	/	DR01B2L-801	CSFR144RSTAY167K4C10


CU7C / CU77C | 6,35 x 3,175 x 2,78 mm flaged  
Applications

HIGH SPEED BEARINGS

Handpiece brand	Handpiece model	Front Bearing (burs side)	Rear Bearing (push side)
CASTELLINI	CLEAN AIR	•	•
	CLEAN AIR 2000	•	•
DABI-ATLANTE	RS350	•	•
DME	ADEN	•	•
FARO	BTC77	•	•
MICRO MEGA	350	•	•
MIDWEST	QUIET AIR LEVER/PUSH BUTTON	•	
	QUIET AIR STANDARD	•	
	TRADITION L/PB CONVERSION	•	
	TRADITION STANDARD	•	
NSK	XGT/TRADITION PB (OEM)	•	•
	KINETIC VIPER 360 TORQUE	•	•
	PHATELUS MINI	•	•
PHATELUS	PHATELUS STND/KINETIC VIPER 360	•	•
	PHATELUS TORQUE	•	•
W & H (ADEC)	ADEC 5000	•	•

CL01LN | 6,35 x 3,175 x 2,78 mm smooth (integral shold)


HIGH SPEED BEARINGS

Size (mm)	
D	6.35
d	3.175
B	2.78

TECNOMED ITALIA Order code	Brand	Cage	Balls
CL01LN	Tecnomed Italia	torlon	ceramic

CROSS REFERENCE BRAND CODE				
Brand	Tecnomed Italia	Myonic	Timken	Barden
Ceramic bearing code	CU77CC	/	DR01B2L-801	CSFR144RSTAY167K4C10

CL01LN | 6,35 x 3,175 x 2,78 mm smooth (integral shold)  
Applications


HIGH SPEED BEARINGS

Handpiece brand	Handpiece model	Front Bearing (burs side)	Rear Bearing (push side)
FARO	DELTA	•	•
	TRILOGY	•	•
	TRILOGY PLUS	•	•
	S400	•	•
	S405	•	•
W & H (ADEC)	700 SERIES	•	•
	800 SERIES	•	•
	TREND TC-80BC	•	•
	TREND TC-95BC	•	•
	TREND TE-95B	•	•
	SYNEA TA-97 LED	•	•
	SYNEA TA-98 LED (Push with 3 cuts)	•	•
	SYNEA TA-98 (Push with 2 cuts)	•	•
	TA-96 MINI	•	•
	VISION TK 100LM	•	•
	190	•	•
	698-796-798-896-898	•	•
	TC-83, TC-98	•	•
	TA-96 L/LC/CM/LCM	•	•
	TA-98 LCM	•	•
	TA-97 CLM/LN	•	•
	TG-97 L, TG-98 L	•	•
	TK-97 LM, TK-98 LM, TK-100 LM	•	•
	E679	•	•
	6500BR	•	•
609	•	•	
TE-95	•	•	
TE-96 LB/RM	•	•	

W & H (ADEC)	TE-97 BC/M/CM/CN	•	•
	TE-98 B/M/CB	•	•
MTC	7807	•	•
TKD	MICRA	•	•
	BRAVIA	•	•
	THERA L-R-E	•	•
BIEN AIR	TITAN	•	•
	ONDINE	•	•
	PRESTIGE	•	•
	S30	•	•
CASTELLINI	CROMA (SMALL HEAD)	•	•
	S32	•	•
	HI-POWER 2 CERANICS	•	•
KAVO	TITANIUM GOLD 2L	•	•
	SILENT POWER 2L/4L/GOLD	•	•
	655-660-S619L	•	•
	8000/B	•	•
SIRONA	E680L	•	•
	9000	•	•
	TI-T2 CONTROL (Push with 2 cuts)	•	•
B.A. INTERNATIONAL	T2 RACER	•	•
	SIROBOOST S/T	•	•
	BA670LKS (Push with 2 cuts)	•	•
	BA695L	•	•
FONA	BA678L	•	•
	8080L	•	•

CN01B / CN11BC / CN11T / CB11BC | 6,35 x 3,175 x 2,78 mm smooth

HIGH SPEED BEARINGS

Size (mm)	
D	6.35
d	3.175
B	2.78

TECNOMED ITALIA Order code	Brand	Cage	Balls
CN01B	Tecnomed Italia	torlon	steel
CN11BC	Tecnomed Italia	torlon	ceramic
CN11T	Myonic	phenolic	steel
CB11BC	Myonic	phenolic	ceramic

CROSS REFERENCE BRAND CODE				
Brand	Tecnomed Italia	Myonic	Timken	Barden
Steel bearing code	CN01B	CN11T	DR55B2L	SR144STAY64K4C10CJ29
Ceramic bearing code	CN11BC	CB11BC	DR55B2L-801	CSR144STAY64K4C10 CJ

CN01B / CN11BC / CN11T / CB11BC | 6,35 x 3,175 x 2,78 mm smooth  
Applications


HIGH SPEED BEARINGS

Handpiece brand	Handpiece model	Front Bearing (burs side)	Rear Bearing (push side)
FARO	DELTA	•	•
	TRILOGY	•	•
	TRILOGY PLUS	•	•
	S400	•	•
	S405	•	•
	700 SERIES	•	•
W & H (ADEC)	800 SERIES	•	•
	TREND TC-80BC	•	•
	TREND TC-95BC	•	•
	TREND TE-95B	•	•
	SYNEA TA-97 LED	•	•
	SYNEA TA-98 LED (Push with 3 cuts)	•	•
	SYNEA TA-98 (Push with 2 cuts)	•	•
	TA-96 MINI	•	•
	VISION TK 100LM	•	•
	190	•	•
MTC	698-796-798-896-898	•	•
	7807	•	•
TKD	MICRA	•	•
	BRAVIA	•	•
	THERA L-R-E	•	•
	TITAN	•	•

BIEN AIR	ONDINE	•	•
	PRESTIGE	•	•
	S30	•	•
	CROMA (SMALL HEAD)	•	•
	S32	•	•
CASTELLINI	HI-POWER 2 CERAMICS	•	•
	TITANIUM GOLD 2L	•	•
	SILENT POWER 2L/4L/GOLD	•	•
KAVO	655-660-S619L	•	•
	8000/B (rotor compatible TMI)	•	•
	E680L	•	•
	9000	•	•
SIRONA	TI-T2 CONTROL (Push with 2 cuts)	•	•
	T2 RACER	•	•
	SIROBOOST S/T	•	•
B.A. INTERNATIONAL	BA670LKS (Push with 2 cuts)	•	•
	BA695L	•	•
	BA678L	•	•
	BA755L	•	•
	BA758L	•	•
FONA	8080L	•	•
IVORY-DENT	T501/L-SM502/4M	•	•

CU6K / CU6CC / CUB6B / CUC99K | 6,35 x 3,175 x 2,78 mm groove

HIGH SPEED BEARINGS

Size (mm)	
D	6.35
D'	6.00
d	3.175
B	2.78
B'	0.88

TECNOMED ITALIA Order code	Brand	Cage	Balls
CU6K	Tecnomed Italia	torlon	steel
CU6CC	Tecnomed Italia	torlon	ceramic
CUB6B	Myonic	phenolic	steel
CUC99K	Myonic	phenolic	ceramic

CROSS REFERENCE BRAND CODE				
Brand	Tecnomed Italia	Myonic	Timken	Barden
Steel bearing code	CU6K	CUB6B	DR70B2L	SR144STAY85K3C10CJ29
Ceramic bearing code	CU6CC	CUC99K	DR70B2L-801	CSR144STAY85K3C10GJ2

CU6K / CU6CC / CUB6B / CUC99K | 6,35 x 3,175 x 2,78 mm groove  
Applications

HIGH SPEED BEARINGS


Handpiece brand	Handpiece model	Front Bearing (burs side)	Rear Bearing (push side)
KAVO	625	•	•
	630	•	•
	635	•	•
	636	•	•
	640	•	•
	640 FG	•	•
	646	•	•
	647	•	•
	649	•	•
	650	+ Code HKV214	•
	BELLA TORQUE FG 642	•	•
	BELLA TORQUE PB 642	•	•
	BELLA TORQUE PB 643	•	•
	CONTACT AIR 632 PB/Screw-type	•	•
MAGNO 634	•	•	
CHIRANA	SATUR TG536	•	•
	SATUR TG546	•	•
	SATUR TGL542	•	•
	SATUR TGL546	•	•
MTI	LYNX PUSH BUTTON	•	•

SIRONA	TS1/2	•	•
	4000S	•	•
	T1 CONTROL (Push with 3 cuts)	•	•
	TC3	•	•
MK-DENT	All products except PRIME LINE & ECO LINE	•	•

CU60SK / CU60CC | 6,35 x 3,175 x 2,35/2,78 mm groove

HIGH SPEED BEARINGS

Size (mm)	
D	6.35
D <sup>1</sup>	6.00
d	3.175
B	2.35
B <sup>1</sup>	0.97
B <sup>2</sup>	2.78

TECNOMED ITALIA Order code	Brand	Cage	Balls
CU60SK	Tecnomed Italia	torlon	steel
CU60CC	Tecnomed Italia	torlon	ceramic

CROSS REFERENCE BRAND CODE				
Brand	Tecnomed Italia	Myonic	Timken	Barden
Steel bearing code	CU60SK	/	DR21J2L-801	/
Ceramic bearing code	CU60CC	/	DR21J2L-801	/


CU60SK / CU60CC | 6,35 x 3,175 x 2,35/2,78 mm groove  
Applications

HIGH SPEED BEARINGS

Handpiece brand	Handpiece model	Front Bearing (burs side)	Rear Bearing (push side)
KAVO	637 B/C	•	•
	E680		•
	6000	•	
	635		•
	637		•
	E675	•	
	E677	•	
	4500	•	•
	5000	•	•
SIRONA	T1 MINI	•	•
MK-DENT	PRIME LINE HC8022	•	•
	PRIME LINE HC9022	•	•

CU700K / CU700KC | 6,35 x 3,175 x 2,38/2,78 mm smooth


HIGH SPEED BEARINGS

Size (mm)	
D	6.35
d	3.175
B	2.78
B'	2.38

TECNOMED ITALIA Order code	Brand	Cage	Balls
CU700K	Tecnomed Italia	torlon	steel
CU700KC	Tecnomed Italia	torlon	ceramic

CROSS REFERENCE BRAND CODE				
Brand	Tecnomed Italia	Myonic	Timken	Barden
Steel bearing code	CU700K	/	DR02J2L	SR144RSTAY196K3
Ceramic bearing code	CU700KC	/	DR02J2L-814	CSR144RSTAY196K3


CU700K / CU700KC | 6,35 x 3,175 x 2,38/2,78 mm smooth  
Applications

HIGH SPEED BEARINGS


Handpiece brand	Handpiece model	Front Bearing (burs side)	Rear Bearing (push side)
KAVO	6000M		•
	7000B (with original rotor)	•	•
	6500		•
	E675		•
	E677		•
	S605	•	•
	S615	•	•

CU58S | 6,35 x 3,175 x 2,78 mm flanged

HIGH SPEED BEARINGS

Size (mm)	
D	6.35
d	3.175
D'	7.50
B	2.78
C'	0.40

TECNOMED ITALIA Order code	Brand	Cage	Balls
CU58S	Tecnomed Italia	torlon	steel

CROSS REFERENCE BRAND CODE				
Brand	Tecnomed Italia	Myonic	Timken	Barden
Steel bearing code	CU58S	/	DR54B2G	/

CU58S | 6,35 x 3,175 x 2,78 mm flanged  
Applications

HIGH SPEED BEARINGS

Handpiece brand	Handpiece model	Front Bearing (burs side)	Rear Bearing (push side)
STAR	430	•	•
	ADVANTAGE	•	•
	VISTA	•	•