


Endo tip for scaler E18, compatible Satelec, DTE, NSK

Endo tip for scaler E18, compatible Satelec, DTE, NSK

Available

SKU: PG130

Categories: Endodontic inserts



Description

Endo tip for scaler E18, compatible Satelec, DTE, NSK.

This Endo tip is used to determine the position of the root canal, remove the dentin, calcification and loose crown at the 1/3 crown of the root canal.

Features: The tip has 12 micro-washers, and it is of good cutting performance.

Back and forth motion of tip brings more efficient operation. In the case of waterless operation, select the mode of "Work for 5~10 seconds and then stop for 2 seconds", and open water every 60 seconds for cooling and rinsing.