



## Fi-P handpiece for vertical condensation, guttapercha system Woodpecker

Fi-P handpiece, Woodpecker Gutta percha system, for hot vertical condensation fillings, used in the root canal filling phase in endodontic treatment.

**Available**

**SKU:** GG210P

**Categories:** Gutta-Percha Devices

### Description

---

Woodpecker Fi-P handpiece is ideal for softening the root canal filling material, the vertical diffusion of the root canal filling material and the compaction of the filling material in the tooth canal.

**How it works** – Fi-P is used for the root canal obturation phase, during the endodontic treatment, by applying the vertical compaction technique.

**Woodpecker Fi-P handpiece features:**

- **Wireless Convenience:** Experience complete freedom of movement during procedures with this fully wireless handpiece.
- **Vertical Condensation:** Specifically indicated for vertical condensation, ensuring optimal adaptation of gutta-percha within the canal.
- **High Precision Temperature Control:** Equipped with advanced temperature regulation for consistent, reliable performance.
- **Rapid Heating:** Heats up in just 0.2 seconds, allowing for efficient, uninterrupted workflow.
- **Adjustable Temperatures:** Features four temperature settings to accommodate various gutta-percha types, ensuring versatility in treatment.
- **Precision Tips:** Includes a 0.35 mm pre-curved tip, designed to penetrate even the narrowest root canals, enhancing treatment outcomes.
- **Long Battery Life:** Powered by a high-performance rechargeable battery, offering up to 4 hours of continuous use or approximately 1500 operations on a full charge.
- **Lightweight Design:** Ergonomically designed to reduce hand fatigue, allowing for comfortable handling during extended procedures.
- **Dual Display for Versatility:** Built-in dual display enables easy manipulation for both right- and left-handed users, enhancing usability.



- **Wide Tip Compatibility:** Compatible with 9 different tip models (WP3504, WP4004, WP4504, WP5004, WP5506, WP5508L, WP5508, WP5510, and WP6012) to meet the needs of various clinical scenarios.
- **Universal Gutta-Percha Compatibility:** Can be used with all brands of gutta-percha, as the heat from the tip efficiently melts the gutta-percha inside the canal.

#### Technical Features

- Weight 80 gr.
- Weight of the charging base 195 gr.
- Preset temperatures: 150°C, 180°C, 200°C and 230°C.
- Heating time in 0.2 seconds.
- Power input AC100V–240V, 50/60 Hz 0.8mA.
- Output power DC15V / 1.6A.
- Battery capacity 2000 mAh.
- Charging time 2.5 hours.
- Power adapter AC / DC.