



OKKIO support arm to be applied on Easy Tray

3 joints articulated arm for EASY TRAY application, for OKKIO medical camera.

Shipping expected in approximately 10 days

SKU: DE1003X.6K

Categories: OKKIO Mounts

Description

The OKKIO 3-joint articulated arm is designed for **exclusive installation on Easy Tray electrified medical trolleys with 5 shelves**, providing a stable, flexible and functional solution for positioning OKKIO cameras or other compatible cameras.

The articulated arm system allows the camera to be positioned at an optimal distance of approximately **300–500mm from the oral cavity**, following the natural line of sight of the clinician during dental procedures. This configuration makes it possible to **capture photos and videos easily and efficiently**, without interfering with clinical workflow.

The high degree of movement ensures high-quality images and stable video recording, eliminating shadow areas and reducing the need for additional personnel or procedural interruptions. Dental professionals can therefore document treatments independently while maintaining continuity and precision during clinical activities.

The OKKIO articulated arm for Easy Tray can be configured according to operational requirements and the technical characteristics of the clinical environment. Thanks to a dedicated adapter, it is also compatible with **commercial cameras equipped with a standard photographic mount**, supporting a **maximum weight of 500g**, further extending its versatility when used on the Easy Tray trolley.