



Scaling tip for scaler, compatible KaVo Piezolux 9, Sonosoft 9

Sickle-like tip for removing supragingival tartar in all quadrants. KaVo Piezolux 9, Sonosoft 9 compatible

Available on request

SKU: PG013K

Categories: Professional Equipment, KaVo

Description

Scaling tip for scaler PG013K.

Detartrasis, or the removal of tartar deposits and supragingival, subgingival and interdental stains, is one of the most common operations performed by all dentists and dental hygienists.

Tecnomed Italia inserts are made from extremely resistant materials and are ideal for use at the high power levels necessary to remove stubborn tartar effectively.

Tips designed for subgingival uses are smaller in size and produce smaller oscillations.

This facilitates work in anatomically inaccessible areas and helps ensure the patient's comfort.

Spray is perfectly tangential to the tip to ensure effective cooling of the operating area.

Sickle-like tip for removing supragingival tartar in all quadrants.

Compatibility Scaling tip for scaler PG013K

- KaVo Piezolux 9;
- Sonosoft 9 compatible.