



Scaling tip for scaler, compatible Satelec 2, DTE GD6, NSK G2

Scaling tip for removing extensive areas of tartar and cement remnants. Satelec 2, DTE GD6, NSK G2 compatible

Available

SKU: PG089

Categories: Professional Equipment, Satelec - DTE - NSK

Description

Scaling tip for scaler PG089.

Detartrasis, or the removal of tartar deposits and supragingival, subgingival and interdental stains, is one of the most common operations performed by all dentists and dental hygienists.

Tecnomed Italia inserts are made from extremely resistant materials and are ideal for use at the high power levels necessary to remove stubborn tartar effectively.

Tips designed for subgingival uses are smaller in size and produce smaller oscillations.

This facilitates work in anatomically inaccessible areas and helps ensure the patient's comfort.

Spray is perfectly tangential to the tip to ensure effective cooling of the operating area.

Scaling tip for removing extensive areas of tartar and cement remnants.

Compatibility Scaling tip for scaler PG089

- Satelec 2;
- DTE GD6;
- NSK G2.