



Scaling tip for scaler, KaVo Piezoled 202 compatible

Scaling tip for supragingival and subgingival tartar removal in all quadrants. Compatible KaVo Piezoled 202

Available

SKU: PG022K

Categories: KaVo

Description

Scaling tip for scaler, compatible KaVo Piezoled 202.

Detartrasis, or the removal of tartar deposits and supragingival, subgingival and interdental stains, is one of the most common operations performed by all dentists and dental hygienists.

Tecnomed Italia inserts are made from extremely resistant materials and are ideal for use at the high power levels necessary to remove stubborn tartar effectively.

Tips designed for subgingival uses are smaller in size and produce smaller oscillations.

This facilitates work in anatomically inaccessible areas and helps ensure the patient's comfort.

Spray is perfectly tangential to the tip to ensure effective cooling of the operating area.

Scaling tip for supragingival and subgingival tartar removal in all quadrants.

Compatibility Scaling tip for scaler: KaVo Piezoled 202.