



Scaling tip for scaler, Mectron S1 compatible

Scaling tip for supragingival tartar removal in all quadrants. Mectron S1 compatible

Available on request

SKU: PG068M

Categories: Scaling inserts

Description

Scaling tip for scaler PG068M.

Detartrasis, or the removal of tartar deposits and supragingival, subgingival and interdental stains, is one of the most common operations performed by all dentists and dental hygienists.

Tecnomed Italia inserts are made from extremely resistant materials and are ideal for use at the high power levels necessary to remove stubborn tartar effectively.

Tips designed for subgingival uses are smaller in size and produce smaller oscillations.

This facilitates work in anatomically inaccessible areas and helps ensure the patient's comfort.

Spray is perfectly tangential to the tip to ensure effective cooling of the operating area.