



Scaling tip for scaler, Mectron S3 compatible

Scaling tip for removing extensive areas of tartar and cement remnants. Mectron S3 compatible

Available

SKU: PG070M

Categories: Mectron - Woodpecker, Professional Equipment

Description

Scaling tip for scaler PG070M.

Detartrasis, or the removal of tartar deposits and supragingival, subgingival and interdental stains, is one of the most common operations performed by all dentists and dental hygienists.

Tecnomed Italia inserts are made from extremely resistant materials and are ideal for use at the high power levels necessary to remove stubborn tartar effectively.

Tips designed for subgingival uses are smaller in size and produce smaller oscillations.

This facilitates work in anatomically inaccessible areas and helps ensure the patient's comfort.

Spray is perfectly tangential to the tip to ensure effective cooling of the operating area.

Scaling tip for removing extensive areas of tartar and cement remnants.

Compatibility Scaling tip for scaler PG070M:

Mectron S3