



Comparative data sheet for models with dryer

Code	150TOVEN	240TOVEN	360TOVEN
Dental units served	1+1	2	3
Air flow	95 liters/minutes	140 liters/minutes	195 liters/minutes
Aspirated volume	152 liters	230 liters	304 liters
Tank capacity	38 liters	38 liters	60 liters
Oil-less	●	●	●
Air dryer	●	●	●
Own cabinet			
Twin cylinder stages	1	1	2
Power	750 Watt	1100 Watt	1500 Watt
Noise level	53 db	53 db	53 db
Power supply	220/240V-50/60 Hz	220/240V-50/60 Hz	220/240V-50/60 Hz
Weight	37 kg	43 kg	58 kg
Dimensions (w x l x h)	430 x 400 x 660 mm	450 x 580 x 750 mm	715 x 387 x 735 mm
Warranty	24 months	24 months	24 months
Certification	CE	CE	CE