



## Comparative data sheet for models without dryer

Code	100TOV	150TOV	150TCV	240TOV	360TOV
Dental units served	1	1+1	1+1	2	3
Air flow	70 liters/minutes	95 liters/minutes	95 liters/minutes	140 liters/minutes	195 liters/minutes
Aspirated volume	115 liters	152 liters	152 liters	230 liters	304 liters
Tank capacity	25 liters	38 liters	38 liters	38 liters	60 liters
Oil-less	●	●	●	●	●
Air dryer					
Own cabinet			●		
Twin cylinder stages	1	1	1	1	2
Power	550 Watt	750 Watt	750 Watt	1100 Watt	1500 Watt
Noise level	52 db	53 db	49 db	52 db	52 db
Power supply	220/240V-50/60 Hz	220/240V-50/60 Hz	220/240V-50/60 Hz	220/240V-50/60 Hz	220/240V-50/60 Hz
Weight	24 kg	28 kg	61 kg	34 kg	51 kg
Dimensions (w x l x h)	400 x 400 x 530 mm	400 x 400 x 680 mm	595 x 570 x 800 mm	450 x 450 x 750 mm	715 x 367 x 735 mm
Warranty	24 months	24 months	24 months	24 months	24 months
Certification	CE	CE	CE	CE	CE