

High speed handpiece controllers



Benchtop controls are best for accurate adjustment of air pressure, turbine speed and spray rate (with spray type handpieces). Speed is best controlled by means of a progressive action foot control.

Our benchtop control units are designed to work with all types of high speed handpiece currently on the market.



Air One is a control unit for sprayless high speed handpieces, complete with filter, pressure regulator and foot control.



Air Two is a control unit for high speed handpieces with an integrated spray system, complete with a 250 ml solution tank, filter, spray on/off selector and foot control. Air Two also incorporates controls for adjusting air and spray pressure.

Comparative data sheet

	Product code AR01P	Product code AR01PM	Product code AR01PHC	Product code AR22P	Product code AR22PM	Product code AR22PHC
Type	Benchtop	Benchtop	Benchtop	Benchtop	Benchtop	Benchtop
Spray				●	●	●
Handpiece connector	2/3 way screw type	4 way screw type	4x6mm tube	2/3 way screw type	4 way screw type	4x6/2x4 mm tube
Progressive foot control	0-3 bar	0-3 bar	0-3 bar	0-3 bar	0-3 bar	0-3 bar
Ai supply	4/6 bar	4/6 bar	4/6 bar	4/6 bar	4/6 bar	4/6 bar
Weight (including foot control)	2.9 kg	2.9 kg	2.9 kg	4.3 kg	4.3 kg	4.3 kg
Dimensions (w x l x h)	150x145x225mm	150x145x225mm	150x145x225mm	250x145x220mm	250x145x220mm	250x145x220mm
Warranty	36 months	36 months	36 months	36 months	36 months	36 months