



Straight handpieces

Product code	SL211C	SL211L
Indications	Dental treatments, restorations	Dental treatments, restorations
Ratio	1:1	1:1
Light intensity up to (lux)	/	40,000
Type of spray	Internal	Internal
Type of coupling	Intramatic Universal	Intramatic Universal
Rotation speed (rpm) <small>(micromotor speed set to 40,000 rpm)</small>	40,000	40,000
Handpiece weight (g)	93	95
Bur diameter (mm)	2.334 / 2.350	2.334 / 2.350
Bur length (mm)	21.95 / 44.50	21.95/44.50
Type of lubrication	Oil	Oil
Autoclave resistant	Yes, 135°C	Yes, 135°C
Packing size lxdxh (mm)	160 x 90 x 43	160 x 90 x 43
Packing weight appr. (g)	250	250
Certification	CE 0476	CE 0476