



Air micromotors

Product code	M714C	M714L
Version	Air driven	Air driven
Compatibility	4 way	6 way
Light source / Voltage (vdc)	No	Led / 3.3-3.8
Type of spray	Internal	Internal
Maximum rotation speed (rpm)	20,000	20,000
Power	Air	Air
Speed adjustment	Yes	Yes
Reverse rotation	Yes	Yes
Handpiece rotation	Yes	No
Air consumption (NL/min at 0.28MPa)	50	50
Handpiece weight (g)	110	130
Handpiece diameter (mm)	20	20
Type of lubrication	Oil	Oil
Autoclave resistant handpiece	Yes-135°C	Yes-135°C
Autoclave resistant removable shell	No	No
Packing size lxdxh (mm)	160 x 90 x 43	160 x 90 x 43
Packing weight appr. (g)	320	320
Certification	CE 0476	CE 0476