

Cabinet for children's operating room











Cabinet in 10/10 galvanised steel sheet, RAL 1023. Worktop in scratchproof glass, RAL 1023.

Product code	DE9.21KDX	DE9.21KSX
Washbasin cabinet	●	●
4 compartment drawer unit	●	●
4 compartment service unit	●	●
Full door unit	1 glass shelf	1 glass shelf
Washbasin – material	Right – glass	Left – glass
Lever type medical tap	●	●
Worktop dimensions mm	2080 x 500 mm	2080 x 500 mm
Height from floor	900 mm	900 mm
Warranty	24 months	24 months

Cabinet for adult operating room



Cabinet in 10/10 galvanised steel sheet, finished in RAL 9010. Worktops are available in two versions: scratchproof glass or water-repellent, scratch-proof laminate finished in RAL 9010. Choice of anodised aluminium handles in eight colours.

Colour		Colour	
	A1 - RAL 5002 Pervinca blue		A8 - RAL 3020 Traffic red
	A2 - RAL 5005 Signal blue		A9 - RAL 5003 Sapphire blue
	A3 - RAL 6027 Light green		A10 - RAL L5018 Acid green
	A6 - RAL 1007 Narcissus yellow		A11 - RAL L3012 Beige red



Product code	DE9.21DX	DE9.21SX	DE9.21DXV	DE9.21SXV	DE9.24DX	DE9.24SX	DE9.25DX	DE9.25SX
Top	Laminate	Laminate	Glass	Glass	Laminate	Laminate	Laminate	Laminate
Washbasin cabinet	●	●	●	●	●	●	●	●
4 compartment drawer unit	●	●	●	●	●	●	●	●
4 compartment service unit	●	●	●	●	●	●	●	●
Full door unit	●1 glass shelf	●1 glass shelf	●1 glass shelf	●1 glass shelf	●1 glass shelf	●1 glass shelf	●with refrigerator	●with refrigerator
Washbasin - material	Right - ABS	Left - ABS	Left - glass	Right - glass	Right - ABS	Left - ABS	Right - ABS	Left - ABS
Lever type medical tap	●	●	●	●	●	●	●	●
Worktop dimensions	2080 x 500 mm	2080 x 500 mm	2080 x 500 mm	2080 x 500 mm	2080 x 500 mm	2080 x 500 mm	2080 x 500 mm	2080 x 500 mm
Height from floor	900 mm	900 mm	900 mm	900 mm	900 mm	900 mm	900 mm	900 mm
Electrovalve					●	●	●	●
Warranty	24 months	24 months	24 months	24 months	24 months	24 months	24 months	24 months
CE	●	●	●	●	●	●	●	●