

Standard Line stools




Smile



Silver


Comfortable and ergonomic, our Smile and Silver stools respect your body's anatomy and stand out for their excellent quality-price ratio. The height adjustment levers are easy to reach and let you adjust seat position without having to interrupt treatment. Our range includes eight different colours for a perfect match with the finish and materials of your dental chair.



Soft, anatomic padding ensures maximum comfort even during long operations. Seamless, fire-repellent upholstery makes these stools easy to clean and disinfect.











All our stools are fitted with pneumatic height adjustment. Our Silver stool also features backrest angle adjustment.










A knob at the back of the seat lets you adjust backrest height on the Silver stool.

Colour and product code table

Available colours	Product code Smile	Product code Silver
 A1 - RAL5002 Periwinkle blue	DE8.02.A1	DE8.01.A1
 A2 - RAL5005 Signal Blue	DE8.02.A2	DE8.01.A2
 A9 - RAL5003 Sapphire blue	DE8.02.A9	DE8.01.A9
 A3 - RAL6027 Light green	DE8.02.A3	DE8.01.A3

Available colours	Product code Smile	Product code Silver
 A10 - RAL5018 Acid green	DE8.02.A10	DE8.01.A10
 A8 - RAL3020 Traffic red	DE8.02.A8	DE8.01.A8
 A11 - RAL3012 Beige red	DE8.02.A11	DE8.01.A11
 A6 - RAL1007 Narcissus yellow	DE8.02.A6	DE8.01.A6

Comparative data sheet for stools

Order code*1	ERGO A	EGGR	StingST	Beauty	Round	Smile	Silver
							
	DE8.06TX	DE8.03TX	DE8.05TX	DE8.07X	DE8.08X	DE8.02X	DE8.01X
Base	Aluminium	Polyurethane	Aluminium	Chrome metal sheet	Chrome metal sheet	ABS	Chrome metal sheet
Seat height adjustment	●	●	●	●	●	●	●
Seat angle adjustment			●				
Backrest height adjustment			●		●		●
Backrest angle adjustment			●		●		●
Backrest depth adjustment				●	●		
Backrest rotation	●	●					
Seat dimensions (mm)	Ø 416	Ø 380	400 x 370	450 x 500	Ø 400	Ø 350	440 x 360
Minimum height (mm)	525	470	530	540	520	450	450 *2
Maximum height (mm)	715	650	720	720	700	600	600 *2
Height extension 150 mm							Optional Product code DE8.014*2
Wheel material	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic
Antistatic wheels							
Weight (kg)	22	22	19	10	10	8	10
Maximum permissible load (kg)	135	135	135	150	150	115	115
Warranty (months)	12	12	12	12	12	12	12
Certificazione	CE	CE	CE	CE	CE	CE	CE

*1 Specify upholstery colour on ordering.

*2 Silver stools with the height extension (DE8.014) have a minimum height of 600 mm and a maximum height of 750 mm.