

# Handy Line stools



Beauty



Round

An appropriate stool is an important investment for the comfort and health of the dental professional.

Beauty and Round professional dentists' chairs are designed for complete comfort even after many hours of work. Their ergonomic structure helps you to maintain correct posture while distributing weight in a uniform manner. Their antibacterial and water-repellent imitation leather upholstery makes these chairs extremely easy to disinfect.

The range offers a choice of thirty different colours for a perfect match with the finish and materials of your dental chair.



Top quality imitation leather upholstery ensures a soft, anatomical seating position with uniform weight distribution and no painful compressions.



Beauty and Round chairs feature adjustable pneumatic height adjustment and backrest depth adjustment.



The Round model also features backrest height and angle adjustment.

## Available colours

RAL 3015 Light pink G17	RAL 4006 Traffic purple G18	RAL 4008 Signal violet G19	RAL 4007 Purple violet G20	RAL 4004 Claret violet G11	RAL 4001 Red lilac G59	RAL 1032 Broom yellow G8	RAL 2011 Deep orange G9	RAL 3031 Orient red G10	RAL 6019 Pastel green G13	RAL 6018 Yellow green G14	RAL 6017 May green G83	RAL 6011 Reseda green G15	RAL 6020 Chrome green G16	RAL 8019 Grey brown G47

RAL 1001 Beige G30	RAL 1002 Sand yellow G31	RAL 8002 Signal brown G32	RAL 9003 Signal white G33	RAL 3020 Traffic red G34	RAL 3004 Purple red G35	RAL 9001 Cream G36	RAL 6027 Light green G37	RAL 6000 Patina green G38	RAL 5024 Pastel blue G39	RAL 5023 Distant blue G40	RAL 5003 Saphire blue G41	RAL 7040 Window grey G42	RAL 7011 Iron grey G43	RAL 9004 Signal black G44

## Order code

Order code*	DE8.07X	DE8.08X
Model	Beauty	Round

\*Specify upholstery colour on ordering.

# Comparative data sheet for stools

	<b>ERGO A</b>	<b>EGG R</b>	<b>StingST</b>	<b>Beauty</b>	<b>Round</b>	<b>Smile</b>	<b>Silver</b>
<b>Order code*¹</b>	<b>DE8.06TX</b>	<b>DE8.03TX</b>	<b>DE8.05TX</b>	<b>DE8.07X</b>	<b>DE8.08X</b>	<b>DE8.02X</b>	<b>DE8.01X</b>
Base	Aluminium	Polyurethane	Aluminium	Chrome metal sheet	Chrome metal sheet	ABS	Chrome metal sheet
Seat height adjustment	●	●	●	●	●	●	●
Seat angle adjustment			●				
Backrest height adjustment			●		●		●
Backrest angle adjustment			●		●		●
Backrest depth adjustment				●	●		
Backrest rotation	●	●					
Seat dimensions (mm)	Ø 416	Ø 380	400 x 370	450 x 500	Ø 400	Ø 350	440 x 360
Minimum height (mm)	525	470	530	420	420	450	450 *²
Maximum height (mm)	715	650	720	560	560	600	600 *²
Height extension 150 mm							Optional Product code DE8.014*²
Wheel material	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic
Antistatic wheels							
Weight (kg)	22	22	19	10	10	8	10
Maximum permissible load (kg)	135	135	135	150	150	115	115
Warranty (months)	12	12	12	12	12	12	12
Certificazione	CE	CE	CE	CE	CE	CE	CE

\*¹ Specify upholstery colour on ordering.

\*² Silver stools with the height extension (DE8.014) have a minimum height of 600 mm and a maximum height of 750 mm.