



# Prophylaxis & Implantology contra-angles

Product code	SL411C	SL411L	SL421C	SL421L	SL415C	SL415L	SL4201C
Indications	Prophylaxis	Prophylaxis	Prophylaxis	Prophylaxis	Restorative	Restorative	Implantology
Ratio	1:1	1:1	2:1	2:1	1:5	1:5	20:1
Light intensity up to (lux)	/	40,000	/	40,000	/	40,000	/
Locking mechanism	Button	Button	Button	Button	Button	Button	Button
Type of spray	Internal	Internal	Internal	Internal	Internal	Internal	External (2)
Type of coupling	Intramatic universal	Intramatic universal	Intramatic universal	Intramatic universal	Intramatic universal	Intramatic universal	Intramatic universal
Rotation speed (rpm) (micromotor speed set to 40,000 rpm)	40,000	40,000	20,000	20,000	200,000	200,000	2,000
Torque (Ncm)	/	/	/	/	/	/	≥ 55
Removable head	Yes	Yes	Yes	Yes	No	No	No
Rotatable head	No	No	No	No	No	No	No
Head diameter (mm)	9.5	9.5	9.5	9.5	9.9	9.9	9.5
Head height (mm)	11.6	11.6	11.6	11.6	14.5	14.5	13.5
Handpiece weight (g)	79	79	79	79	82	82	58
Bur/file diameter (mm)	2.334/2.350	2.334/2.350	2.334/2.350	2.334/2.350	1.59/1.60	1.59/1.60	2.334/2.350
Bur/file length (mm)	21.95	21.95	21.95	21.95	10/25	10/25	21.95
Type of lubrication	Oil	Oil	Oil	Oil	Oil	Oil	Oil
Autoclave resistant	Yes, 135°C	Yes, 135°C	Yes, 135°C	Yes, 135°C	Yes, 135°C	Yes, 135°C	Yes, 135°C
Packing size lxdxh (mm)	160 x 90 x 43	160 x 90 x 43	160 x 90 x 43	160 x 90 x 43	160 x 90 x 43	160 x 90 x 43	160 x 90 x 43
Packing weight appr. (g)	250	250	250	250	250	250	250
Certification	CE 0476	CE 0476	CE 0476	CE 0476	CE 0476	CE 0476	CE 0476