

AIR-FLOW HANDPIECES



Prophylaxis plays an essential role in preventing pathologies of the oral cavity, but to be effective, it requires the use of user-friendly, professional instruments. Prophy Silk by Ivory air-driven dental polisher uses sodium bicarbonate powder and is delicate on the enamel. The Prophy Silk polisher is ideal for removing plaque, treating discoloration and for general cleaning. The principle of functioning is based on the mechanical action of a flow of sodium bicarbonate crystals that is accelerated in a jet of compressed air and water and dissolved to form a high-power spray. The jet expelled from the tip produces a powerful cleaning effect and is pleasant for the patient too. This easy to use instrument is ideal for dentists and dental hygienists alike.



Thanks to the quick connector for the cord of your dental unit, you can perform consecutive prophylactic treatments more rapidly and achieve excellent results on all your patients. Separate instruments with trailing power supplies will be a distant memory.

Unrivalled lightness

Reduced weight, compact design and an ergonomic narrowing in the body give the handpiece an excellent natural grip and help you to adopt an ergonomically correct posture, avoiding the hand cramps caused by unnatural hand positions. The hermetically sealed sodium bicarbonate container is made from a tough polymer and can be refilled with powder easily and safely. The handpiece can be disinfected in an autoclave at up to 135°C for maximum hygiene.

Maximum hygiene

A special nozzle design generates a jet in the form of an inverted champagne glass. This avoids powder dispersal and ensures continuous rinsing of the oral cavity. The tip can rotate through 360° to permit rapid and effective cleaning even in areas of difficult access. The tip can be disinfected in an autoclave at up to 135°C for maximum hygiene.



**Handpiece
REF. T501LK**

**Handpiece
REF. T501LN**

**Handpiece
REF. T501LS**

**Handpiece
REF. T501LB**

**Handpiece
REF. S502B**