








Air-flow handpieces

Product code	P900B	P900K	P900S	P900N	P900W
Compatibility	Bien Air®	Kavo®	Sirona®	NSK®	W6H®
Instrument connection					
360° tip rotation	Yes	Yes	Yes	Yes	Yes
Handpiece weight (g)	356	356	356	356	356
Autoclave resistant (only tip)	Yes	Yes	Yes	Yes	Yes
Packing size lxdxh (mm)	190 x 125 x 75	190 x 125 x 75	190 x 125 x 75	190 x 125 x 75	190 x 125 x 75
Packing weight appr. (g)	500	500	500	500	500
Certification	CE	CE	CE	CE	CE