



Silent Line compressors



The new Silent Line range of compressors features a modern design that makes them extremely silent running. The new version features a latest generation air dryer that is superbly efficient and guarantees a perfectly dry air supply to your dental unit's air lines. Because Tecnomed Italia compressors deliver top quality, high purity compressed air, they help ensure the correct functioning of instruments and reduce maintenance costs.



Comparative data sheet for models with dryer

Code	150TOVEN	240TOVEN	360TOVEN
Dental units served	1+1	2	3
Air flow	95 liters/minutes	140 liters/minutes	195 liters/minutes
Aspirated volume	152 liters	230 liters	304 liters
Tank capacity	38 liters	38 liters	60 liters
Oil-less	●	●	●
Air dryer	●	●	●
Own cabinet			
Twin cylinder stages	1	1	2
Power	750 Watt	1100 Watt	1500 Watt
Noise level	53 db	53 db	53 db
Power supply	220/240V-50/60 Hz	220/240V-50/60 Hz	220/240V-50/60 Hz
Weight	37 kg	43 kg	58 kg
Dimensions (w x l x h)	430 x 400 x 660 mm	450 x 580 x 750 mm	715 x 387 x 735 mm
Warranty	24 months	24 months	24 months
Certification	CE	CE	CE



Comparative data sheet for models without dryer

Code	100TOV	150TOV	150TCV	240TOV	360TOV
Dental units served	1	1+1	1+1	2	3
Air flow	70 liters/minutes	95 liters/minutes	95 liters/minutes	140 liters/minutes	195 liters/minutes
Aspirated volume	115 liters	152 liters	152 liters	230 liters	304 liters
Tank capacity	25 liters	38 liters	38 liters	38 liters	60 liters
Oil-less	●	●	●	●	●
Air dryer					
Own cabinet			●		
Twin cylinder stages	1	1	1	1	2
Power	550 Watt	750 Watt	750 Watt	1100 Watt	1500 Watt
Noise level	52 db	53 db	49 db	52 db	52 db
Power supply	220/240V-50/60 Hz	220/240V-50/60 Hz	220/240V-50/60 Hz	220/240V-50/60 Hz	220/240V-50/60 Hz
Weight	24 kg	28 kg	61 kg	34 kg	51 kg
Dimensions (w x l x h)	400 x 400 x 530 mm	400 x 400 x 680 mm	595 x 570 x 800 mm	450 x 450 x 750 mm	715 x 367 x 735 mm
Warranty	24 months	24 months	24 months	24 months	24 months
Certification	CE	CE	CE	CE	CE