











# Couplings

Product code	C31K	C30LK	C31LK	C30LW	C30LN	C30S	C30LS	C30LB
Compatibility	Kavo® Multiflex®	Kavo® Multiflex®	Kavo® Multiflex®	W&H® Roto Quick® RQ 24	NSK® PTL-CL-LED	Sirona® Click&Go	Sirona® Click&Go	Bien Air® Unifix
Instrument connection								
Light intensity (lux)	No	25,000	25,000	No	25,000	No	25,000	No
Light source / Voltage (vdc)	No	Led / 3.3-3.8	Led / 3.3-3.8	Yes	Led / 3.3-3.8	No	Led / 3.3-3.8	Yes
360° instrument rotation	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Spray adjustment	Yes	No	Yes	No	No	No	Yes	Yes
Couplings weight (g)	31	24	31	31	31	30	37	30
Autoclave resistant	No	No	No	No	No	No	No	No
Packing size lxdxh (mm)	105 x 75 x 25	105 x 75 x 25	105 x 75 x 25	105 x 75 x 25	105 x 75 x 25	105 x 75 x 25	105 x 75 x 25	105 x 75 x 25
Packing weight appr. (g)	120	120	120	120	120	120	120	120
Certification	CE	CE	CE	CE	CE	CE	CE	CE