

LUBRICANT OIL

The Ivory range of lubricants is universally recognised for its unrivalled degreasing and cleaning power. The elimination of overheating in moving parts, the prevention of deposits, lower consumption and running costs, extended instrument life, zero thermal stress and oxidation, better performance and fewer technical assistance callouts are just some of the benefits of using these oils.

FAST-OIL, THE NEW FRONT LUBRICATION SYSTEM

The innovative front lubrication system **Fast-Oil** designed and developed by Ivory for Silk Line instruments represents one of the most important developments on the market today. Thanks to this clever design, you can spray-lubricate the turbines of Silk Line instruments from the front, without any need for adapters. This is possible thanks to a revolutionary lubricant distribution system inside the rotor. The result is significantly reduced repair costs, less chance of mistakes when servicing the instrument and reduced maintenance times too.



Adapter for turbines
and micromotors
REF. CL5.001V



Adapter for contra-angles
and heads
REF. LIV1013



Silk Oil lubricant oil
REF. OL101S