


Monitor arms



Though they are indispensable items of equipment in any modern dental studio, monitors and computers must not compromise the high levels of hygiene needed. Medical PCs and monitors must therefore conform to strict hygiene standards as well as providing a high level of functionality in terms of mobility and usability. Stylish, attractive looks are also a plus. Our arms support all types of monitor. Installation is extremely easy and does not require specialist assistance, just a simple Allen key! They also have a convenient handle so that you can move the monitor or all in one PC around easily. Fixing holes are provided for Vesa 75 and Vesa 100 mountings. This new range of modular monitor arms from Tecnomed Italia offers amazing versatility and lets you specify components of various sizes. You can therefore configure your monitor arm exactly as you want it. Cables can be routed through the arm as an option. Various cables are available, including VGA, USB and power.


Modular jointed horizontal arm

Yoke




	Code GA459010	Code GA509010	Code GA559010
Colour	Dental white*	Dental white*	Dental white*
L	115 mm	115 mm	115 mm
H	90 mm	90 mm	90 mm
A	ø 35 mm	ø 35 mm	ø 35 mm
B	ø 45 mm	ø 50 mm	ø 55 mm
Hollow yoke	●	●	●
Warranty	24 months	24 months	24 months

Adapter



	Code GA459010	Code GA509010	Code GA559010
Colour	Dental white*	Dental white*	Dental white*
L	115 mm	115 mm	115 mm
H	90 mm	90 mm	90 mm
A	ø 35 mm	ø 35 mm	ø 35 mm
B	ø 45 mm	ø 50 mm	ø 55 mm
Hollow yoke	●	●	●
Warranty	24 months	24 months	24 months


Arm



	Code BR309010	Code BR409010
Colour	Dental white*	Dental white*
L	300 mm	400 mm
H	220 mm	220 mm
A	ø 35 mm	ø 35 mm
Hollow yoke	●	●
Warranty	24 months	24 months

* Colour tones between RAL9003 e RAL9010

VESA bracket



	Code SV019010	Code SV029010
Colour	Dental white*	Dental white*
Handle	●	
A	115 mm	115 mm
B	115 mm	115 mm
L	245 mm	
Warranty	24 months	24 months

Bracket



	Code SV039003
L	363 mm
A	225 mm
B	43 mm
Warranty	24 months

Complete Eco arm



	Code BR45EC
Colour	Dental white*
L	330 mm
Yoke Ø	42-45 mm
Hollow arm	Cable duct
Hollow yoke	NO
VESA bracket	75-100
Warranty	24 months

* Colour tones between RAL9003 e RAL9010

Accessories



5 metre amplified USB cable.
Product code CL04USB



5 metre VGA cable. Connector fitting.
Product code CL05VGA