

Comfort Line stools



ERGO A



EGG R



StingST

An appropriate stool is an important investment for the comfort and health of the dental professional. Tecnomed has a vast experience in the design and production of ergonomic dental stools. Tecnomed stools help you maintain a correct posture with natural curvature of the spine, uniform distribution of body weight and no painful compressions. EGG R, ERGO A and StingST stools are designed to offer superior comfort for dentists and patients at all times.



Ergo A, EGG R and StingST stools are upholstered in top quality imitation leather and are free from stitching for maximum hygiene and easy cleaning. The upholstery material is fireproof according to ISO 6875 and UNI 4818.74. This top quality upholstery helps ensure a soft, anatomical seating position with uniform weight distribution and no painful compressions.



The StingST stool is specially designed for dentists who work in a seated position for many hours at a time. Its special saddle shape seat promotes correct posture. The "V" leg position also relieves strain on the spine and improves circulation. Seat angle is adjustable to suit requirements.



The "U" shaped, rotating backrest on the Ergo A and EGG R models provides good lumbar support and a practical rest for the user's arm. The backrest can be reversed, right to left, making these stools the perfect choice even for left handed dentists.

Colour and product code table

Available colours	Product code ERGO A	Product code EGG R	Product code StingST
 RAL 9017 Black	DE8.06T1007	DE8.03T1007	DE8.05T1007
 RAL 9002 White	DE8.06T9002	DE8.03T9002	DE8.05T9002
 RAL 5003 Dark blue	DE8.06T5003	DE8.03T5003	DE8.05T5003
 RAL 5002 Electric blue	DE8.06T1009D	DE8.03T1009D	DE8.05T1009D
 RAL 5012 Periwinkle blue	DE8.06T3491	DE8.03T3491	DE8.05T3491
 RAL 6000 Leaf green	DE8.06T5511	DE8.03T5511	DE8.05T5511
 RAL 6027 Acqua marine green	DE8.06T3493	DE8.03T3493	DE8.05T3493

Available colours	Product code ERGO A	Product code EGG R	Product code StingST
 RAL 6019 Apple green	DE8.06T5515	DE8.03T5515	DE8.05T5515
 RAL 3003 Red	DE8.06T1414	DE8.03T1414	DE8.05T1414
 RAL 4008 Purple	DE8.06T5510	DE8.03T5510	DE8.05T5510
 RAL 4009 Cyclamen violet	DE8.06T3492	DE8.03T3492	DE8.05T3492
 RAL 3022 Salmon	DE8.06T7791	DE8.03T7791	DE8.05T7791
 RAL 1017 Saffron yellow	DE8.06T5503	DE8.03T5503	DE8.05T5503

Handy Line stools



Beauty



Round

An appropriate stool is an important investment for the comfort and health of the dental professional. Beauty and Round professional dentists' chairs are designed for complete comfort even after many hours of work. Their ergonomic structure helps you to maintain correct posture while distributing weight in a uniform manner. Their antibacterial and water-repellent imitation leather upholstery makes these chairs extremely easy to disinfect. The range offers a choice of thirty different colours for a perfect match with the finish and materials of your dental chair.



Top quality imitation leather upholstery ensures a soft, anatomical seating position with uniform weight distribution and no painful compressions.



Beauty and Round chairs feature adjustable pneumatic height adjustment and backrest depth adjustment.



The Round model also features backrest height and angle adjustment.

Available colours

RAL 3015 Light pink G17	RAL 4006 Traffic purple G18	RAL 4008 Signal violet G19	RAL 4007 Purple violet G20	RAL 4004 Claret violet G11	RAL 4001 Red lilac G59	RAL 1032 Broom yellow G8	RAL 2011 Deep orange G9	RAL 3031 Orient red G10	RAL 6019 Pastel green G13	RAL 6018 Yellow green G14	RAL 6017 May green G83	RAL 6011 Reseda green G15	RAL 6020 Chrome green G16	RAL 8019 Grey brown G47
RAL 1001 Beige G30	RAL 1002 Sand yellow G31	RAL 8002 Signal brown G32	RAL 9003 Signal white G33	RAL 3020 Traffic red G34	RAL 3004 Purple red G35	RAL 9001 Cream G36	RAL 6027 Light green G37	RAL 6000 Patina green G38	RAL 5024 Pastel blue G39	RAL 5023 Distant blue G40	RAL 5003 Sapphire blue G41	RAL 7040 Window grey G42	RAL 7011 Iron grey G43	RAL 9004 Signal black G44

Order code

Order code*	DE8.07X	DE8.08X
Model	Beauty	Round

*Specify upholstery colour on ordering.

Standard Line stools




Smile




Silver


Comfortable and ergonomic, our Smile and Silver stools respect your body's anatomy and stand out for their excellent quality-price ratio. The height adjustment levers are easy to reach and let you adjust seat position without having to interrupt treatment. Our range includes eight different colours for a perfect match with the finish and materials of your dental chair.



Soft, anatomic padding ensures maximum comfort even during long operations. Seamless, fire-repellent upholstery makes these stools easy to clean and disinfect.











All our stools are fitted with pneumatic height adjustment. Our Silver stool also features backrest angle adjustment.










A knob at the back of the seat lets you adjust backrest height on the Silver stool.

Colour and product code table

Available colours	Product code Smile	Product code Silver
 A1 - RAL5002 Periwinkle blue	DE8.02.A1	DE8.01.A1
 A2 - RAL5005 Signal Blue	DE8.02.A2	DE8.01.A2
 A9 - RAL5003 Sapphire blue	DE8.02.A9	DE8.01.A9
 A3 - RAL6027 Light green	DE8.02.A3	DE8.01.A3

Available colours	Product code Smile	Product code Silver
 A10 - RAL5018 Acid green	DE8.02.A10	DE8.01.A10
 A8 - RAL3020 Traffic red	DE8.02.A8	DE8.01.A8
 A11 - RAL3012 Beige red	DE8.02.A11	DE8.01.A11
 A6 - RAL1007 Narcissus yellow	DE8.02.A6	DE8.01.A6

Comparative data sheet for stools

Order code*1	ERGO A	EGG R	StingST	Beauty	Round	Smile	Silver
	 DE8.06TX	 DE8.03TX	 DE8.05TX	 DE8.07X	 DE8.08X	 DE8.02X	 DE8.01X
Base	Aluminium	Polyurethane	Aluminium	Chrome metal sheet	Chrome metal sheet	ABS	Chrome metal sheet
Seat height adjustment	●	●	●	●	●	●	●
Seat angle adjustment			●				
Backrest height adjustment			●		●		●
Backrest angle adjustment			●		●		●
Backrest depth adjustment				●	●		
Backrest rotation	●	●					
Seat dimensions (mm)	Ø 416	Ø 380	400 x 370	450 x 500	Ø 400	Ø 350	440 x 360
Minimum height (mm)	525	470	530	420	420	450	450 *2
Maximum height (mm)	715	650	720	560	560	600	600 *2
Height extension 150 mm							Optional Product code DE8.014*2
Wheel material	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic	Plastic
Antistatic wheels							
Weight (kg)	22	22	19	10	10	8	10
Maximum permissible load (kg)	135	135	135	150	150	115	115
Warranty (months)	12	12	12	12	12	12	12
Certificazione	CE	CE	CE	CE	CE	CE	CE

*1 Specify upholstery colour on ordering.

*2 Silver stools with the height extension (DE8.014) have a minimum height of 600 mm and a maximum height of 750 mm.