

# Neo Line

## The value of simplicity

CONTRA-ANGLE AND STRAIGHT  
HANDPIECES



**IVORU**  
Mechanikkonzept

*European reliability  
for your professionalism*  
MADE IN EU

## INNOVATION AND QUALITY FOR ALL YOUR NEEDS

---

Ivory Mechanikkonzept was founded in 2015 by a group of professionals with more than 25 years of experience in the dental sector, quickly establishing itself in more than 15 countries as Smart Company, a manufacturer of high quality dental handpieces "Made in Europe".

With the mission "**Make the doctor's work easier and more effective**", we have activated innovative and exclusive projects such as online training, we use new technologies and cutting-edge manufacturing techniques, we offer the competitive advantage of handpieces compatible with the most famous brands in the market.

## NEO LINE, THE NEW FRONTIER OF EVOLUTION

---

The new Neo Line handpieces range, combines the excellent Ivory quality with a new design, guaranteeing the physician maximum precision and safety at an even more affordable price.



# CONTRA-ANGLE HANDPIECE

The excellence you were looking for is now at your fingertips.

- + **Removable head**  
Quick cleaning and maintenance.
- + **Small size**  
Also reach areas that are difficult to access.
- + **Extreme lightness**  
Work without stress.
- + **Reinforced bearings**  
Excellent and durable performance.
- + **Perfect balance**  
Any intervention in total absence of vibrations.
- + **Internal spray**
- + **Transmission 1:1**  
For any prophylaxis intervention.
- + **Cost optimisation**  
Indeed, it is possible to use different heads, normal or reducers, on the same angle.
- + **Italian design**  
The path of design always goes through Italy.
- + **2 years warranty**  
Quality production
- + **CE0476**  
European medical certificate 0476.

Contra-angle REF. SL411E



LIGHTNESS WITHOUT  
VIBRATIONS



COMPATIBILITY  
WITH  
KAVO® INTRA HEADS



SMALL SIZE HEAD



WEAR RESISTANT  
FRONT  
BEARING

# STRAIGHT HANDPIECE

Made to personalise your work in a highly professional and safe way.

- + **Small size**  
Also reach areas that are difficult to access.
- + **Extreme lightness**  
Work without stress.
- + **Perfect balance**  
Any intervention in total absence of vibrations.
- + **Internal spray**
- + **Transmission 1:1**  
For any prophylaxis intervention.
- + **Italian design**  
The path of design always goes through Italy.
- + **2 years warranty**  
Quality production
- + **CE0476**  
European medical certificate 0476.

Straight handpiece REF. SL211E



MAXIMUM SILENCE



ABSENCE  
OF VIBRATIONS

# TECHNICAL DATA SHEET



Model	SL411E	SL211E
Removable head	YES	-
Optical fibre	NO	NO
Transmission	1:1	1:1
Indications	Prophylaxis	Prophylaxis
Stop mechanism	Button/Push	-
Type of spray	Internal	Internal
Compatibility	Intramatic Universal	Intramatic Universal
Maximum revolutions per minute (RPM)	40,000	40,000
Head diameter	9.5 mm	-
Head height	11.6 mm	-
Cutter diameter	2.334/2.350 mm	2344/2350 mm
Cutter length	21.95 mm	21.95/44.50 mm
Weight	79 gr	93 gr
Autoclavable	YES, 135°	YES, 135°
Type of lubrication	Oil	Oil

Key:- = not available



**IVORY SRL**

**Registered office**

Via Larga, 7  
20122 Milan (MI) - Italy  
[ivory@ivory-dent.com](mailto:ivory@ivory-dent.com)  
[www.ivory-dent.com](http://www.ivory-dent.com)

**IVORY AUSTRIA**

**Zentraler Verkaufsabteilung**

Twin Tower, Wienerbergstrasse 11  
1100 Wien - Österreich  
[austria@ivory-dent.com](mailto:austria@ivory-dent.com)

**IVORY ITALY**

**Sales Office**

via Roma, 9  
61040 Monte Porzio (PU) - Italy  
[italy@ivory-dent.com](mailto:italy@ivory-dent.com)

[www.ivory-dent.com](http://www.ivory-dent.com)