



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/11mm-suction-cannula-adapter-for-skudo-including-o-ing-in-derling/>

11mm suction cannula adapter for SKUDO, including o-ing, in derling

11mm suction cannula adapter for SKUDO, including o-ing, in derling.

11,32 €

Available

SKU: DP2003.11

Categories: Skudo

Description

11mm cannula adapter for SKUDO.

11mm suction cannula adapter for SKUDO, including o-ing, in derling.