



12-piece compatible EMS insert kit, NO branded

12 pieces **EMS compatible inserts kit** for scaling, periodontics and endodontics, with E, P and G series tips for everyday clinical use.

Available

SKU: OFOUT107

Categories: Outlet

Description

12-piece EMS compatible inserts kit consisting of a selection of inserts for ultrasonic scaling, periodontal applications and endodontic procedures.

The set includes E, P and G series tips for versatile and practical clinical use.

Kit composition

Code	Quantity	Description
PG054	1	E3 tip
PG066	3	E1 tip + 6-file kit
PG062	2	P3 tip
PG053	2	P4D tip
PG031	3	G4 tip
PG042	1	G6 tip