

Articulated monitor arm, 300 mm length, AESTHETIC IMPERFECTION



Monitor support arm for dental unit application, with joint and 35mm spigot, 300mm length. **PRODUCT WITH AESTHETIC IMPERFECTIONS ONLY 1 PIECE AVAILABLE.**

Available

SKU: OFOUT107

Categories: Outlet

Description

BR309010: articulated monitor support arm with joint and 35mm spigot, 300mm length. **PRODUCT WITH AESTHETIC IMPERFECTIONS ONLY 1 PIECE AVAILABLE.**

Main features

Vertical pole mounting
Internal routing for electrical cables

Technical specifications

Arm spigot diameter ?35mm
Dimensions L x H: 300 x 190mm
Material: metal
Colour: dental white
Certification: 2001/95/EC
Warranty: 36 months

Application areas

Dental / podiatric treatment units
Equipped poles

Configuration note

The monitor support arm can be combined with VESA mounts, blocks or jaws compatible with a 35mm spigot.