



Articulated monitor arm, 400 mm length, AESTHETIC IMPERFECTION

Articulated monitor support arm with joint and 35mm spigot, 400mm length, **PRODUCT WITH AESTHETIC IMPERFECTION, ONLY 2 PIECES AVAILABLES.**

Available

SKU: OFOUT106

Categories: Outlet

Description

BR409010: articulated monitor support arm with joint and 35mm spigot, 400mm length. – **PRODUCT WITH AESTHETIC IMPERFECTION, ONLY 2 PIECES AVAILABLES.**

Main features

Vertical pole mounting
Internal routing for electrical cables

Technical specifications

Arm spigot diameter ?35mm
Dimensions L x H: 400 x 190mm
Material: metal
Colour: dental white
Certification: 2001/95/EC
Warranty: 36 months

Application areas

Dental / podiatric treatment units
Equipped poles

Configuration note

The monitor support arm can be combined with VESA mounts, blocks or jaws compatible with a 35mm spigot.