



## 13-piece compatible Satelec insert kit, NO branded

**Satelec DTE compatible inserts kit** with 13 ultrasonic tips for scaling and endodontics, ideal for canal treatment, calculus removal and precision clinical procedures.

**Available**

**SKU:** OFOUT106

**Categories:** Outlet

### Description

#### 13 piece Satelec DTE compatible inserts kit for scaling and endodontics

**Satelec compatible inserts kit** designed for everyday ultrasonic scaling and endodontic procedures. This 13-piece set offers a practical combination of compatible tips selected for professional clinical use. The kit includes a balanced assortment for different applications, from routine scaling to more specific endodontic work. It combines standard working tips with dedicated shapes for long canals and 90° access.

#### Kit composition 13 piece Satelec DTE compatible inserts

Code	Qty	Description
PG089	2	G2 tip
PG118	1	E4 tip
PG119	1	E5 tip
PG121	1	Long canals tip
PG122	1	90° endo tip
PG123	1	E10 tip
PG125	1	ED11 tip
PG126	1	E11D tip
PG127	1	E14 tip
PG128	1	E15 tip