



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/3-shelf-surgical-trolley-1000x500x935-mm/>



3-shelf service cart 1010x590x940 mm, shelves 1000x500 mm

Large service cart for dental / medical use with wheels, 3 shelf, made entirely of stainless steel.

421,26 €

Available on request

SKU: DE1005X

Categories: Multifunctional cart

Description

3 shelf surgical trolley – Large service cart for dental / medical use with wheels, 3 shelf, made entirely of stainless steel.

Dimensions (H.L.W) 101x94x59cm. Shelf size: 100x50 cm.