



### 3-shelf service cart 1010x590x940 mm, shelves 1000x500 mm

Large service cart for dental / medical use with wheels, 3 shelf, made entirely of stainless steel.

**Available on request**

**SKU:** DE1005X

**Categories:** Carts

## Description

---

The **stainless-steel medical service trolley** with **three shelves** is designed to provide **efficient organisation and hygienic handling** of instruments and materials in **dental clinics, hospitals and operating-theatre environments**.

Its **fully welded stainless-steel frame** ensures excellent rigidity, durability and easy cleaning, while **smooth, joint-free surfaces** help maintain high hygiene standards.

The **round-tube side supports** give the trolley a solid structure and a refined, functional design suited to modern clinical interiors.

With overall dimensions of **1010 × 590 × 940 mm** and shelves measuring **1000 × 500 mm**, it provides ample working space for trays, tools or accessories.

Each shelf is fitted with a **raised edge** to prevent items from slipping during transport, guaranteeing safe and stable use.

The trolley runs on **swivelling castors** that ensure **quiet and effortless mobility**, allowing smooth movement even in confined areas such as treatment rooms or sterilisation zones.

Compact yet robust, it supports everyday operations where **hygiene, ergonomics and reliability** are essential.

Thanks to its **minimalist design, corrosion-resistant construction and easy-clean finish**, this stainless-steel medical service trolley is ideal for intensive use in **surgical theatres, dental surgeries and clinical laboratories**.

It represents a dependable, professional-grade piece of equipment that enhances workflow efficiency and upholds stringent cleanliness requirements.