



Adjustable make-up chair, backrest tilt 45 degrees, white

Adjustable make-up chair h 55/72 cm, back tilt 45 degrees, sky colour white

Available on request

SKU: ES1005BP

Categories: Beauty & Spa Treatment Furniture

Description

The adjustable white make-up chair combines elegance, comfort and professional functionality, making it ideal for beauty salons, make-up studios and spa environments.

Its clean and elegant design fits perfectly into modern professional settings.

Equipped with a **reclining backrest (up to 45°)** and a **telescopic headrest**, it ensures maximum comfort for the client during make-up or facial treatments.

The **removable armrests and footrest** allow full flexibility for the operator and enhanced comfort for the client.

The seat **rotates 360°** and the **hydraulic version includes a locking mechanism** for maximum safety and precision.

Height adjustment is possible via hydraulic or gas pump, adapting to different working positions.

Built with a **strong metal structure** and **washable eco-leather upholstery**, it is durable, hygienic and easy to maintain.

With a lifting capacity of 150 kg and reclining angle from 0° to 20°, the adjustable white make-up chair offers the perfect balance of stability and comfort.

Optional **padded armrests with manicure bowls** can be added for increased functionality in professional use.