



Printed on 02/06/2025

<https://shop.tecnomeditalia.com/en/product/adjustable-make-up-chair-backrest-tilt-45-degrees-white/>



Adjustable make-up chair, backrest tilt 45 degrees, white

Adjustable make-up chair h 55/72 cm, back tilt 45 degrees, sky colour white

1.579,42 €

Available on request

SKU: ES1005BP

Categories: Armchairs and massage tables

Description

Adjustable make-up chair.

Practical make-up chair with a sober, elegant design.

It is equipped with telescopic head rest, a reclining back (20°), with removable arms and foot rest.

The seat rotates through 360° and has a block on the hydraulic version. The height is adjustable via hydraulic pump action or gas.

The arms are removable and available in two models:

- STANDARD: Black polyurethane.
- OPTIONAL: Padded and fitted with manicure trays.

Technical features Adjustable make-up chair.

Lifting capacity: 150 Kg.

Tilting: 0° – 20°

Dimensions: 45 cm x 100 cm

Height adjustment:

- pump action from 45 cm to 48 cm
- pump gas from 55 cm to 72 cm

Weight: 30 kg