



## Adjustable make-up chair, backrest tilt 45 degrees, white

Adjustable make-up chair h 55/72 cm, back tilt 45 degrees, sky colour white

**Available on request**

**SKU:** ES1005BP

**Categories:** Beauty & Spa Treatment Furniture

### Description

---

The **adjustable white make-up chair** combines **elegance, comfort and professional functionality**, making it ideal for **beauty salons, make-up studios and spa environments**. Its **clean and elegant design** fits perfectly into modern professional settings.

Equipped with a **reclining backrest (up to 45°)** and a **telescopic headrest**, it ensures maximum comfort for the client during make-up or facial treatments.

The **removable armrests and footrest** allow full flexibility for the operator and enhanced comfort for the client.

The seat **rotates 360°** and the **hydraulic version includes a locking mechanism** for maximum safety and precision.

Height adjustment is possible via **hydraulic or gas pump**, adapting to different working positions.

Built with a **strong metal structure** and **washable eco-leather upholstery**, it is durable, hygienic and easy to maintain.

With a **lifting capacity of 150 kg** and **reclining angle from 0° to 20°**, the **adjustable white make-up chair** offers the perfect balance of stability and comfort.

Optional **padded armrests with manicure bowls** can be added for increased functionality in professional use.