



AIRONE Set ? Dental lab turbine control unit with included turbin

Kit including AIRONE 1-way turbine control unit, laboratory turbine and pneumatic foot pedal.

Available on request

SKU: OFR01P

Categories: Professional Equipment

Description

AIRONE Dental Lab Turbine Control Unit with Included Turbine

The **AIRONE turbine control unit** is a bench-top system designed to power and control dental laboratory turbines used for finishing and processing prosthetic materials.

The system allows precise adjustment of **air pressure and turbine speed**, ensuring stable performance during laboratory procedures involving zirconia, dental ceramics, metals and resins.

The **progressive pneumatic foot pedal** enables accurate speed control while keeping both hands free, improving ergonomics and operational precision.

The unit includes an **air pressure regulator and integrated air filter**, ensuring clean and stable airflow while protecting the turbine during operation.

The AIRONE system is supplied with a **laboratory turbine included**, providing a ready-to-use solution for dental laboratories.

Kit includes

- Pneumatic control unit with progressive foot pedal
Tecnomed Italia code **AR01PH**
- Laboratory turbine without spray
Speed **300,000 rpm**
FG bur locking system