



AIRTWO Set ? Dental lab turbine control unit with included turbine

Kit including AIRTWO 4-way turbine control unit, laboratory spray turbine, 250 ml bottle and pneumatic foot pedal for speed control.

Available on request

SKU: OFR22P

Categories: Professional Equipment

Description

AIRTWO Dental Lab Turbine Control Unit with Included Turbine

The **AIRTWO turbine control unit** is a bench-top system designed to power and control spray turbines used in dental laboratories during finishing and processing procedures for prosthetic materials.

The system allows precise adjustment of **air pressure, spray flow and turbine speed**, ensuring stable performance during laboratory procedures involving zirconia, dental ceramics, metals and resins.

The **progressive pneumatic foot pedal** enables accurate turbine speed modulation while keeping both hands free during operation, improving ergonomics and working precision.

The unit includes an **air pressure regulator and air filter**, ensuring clean and stable airflow while helping protect the turbine during use.

The AIRTWO kit is supplied with a **laboratory spray turbine included**, providing a ready-to-use solution for dental laboratories.

Kit includes

- Pneumatic control unit with **250 ml spray bottle** and progressive foot pedal
Tecnomed Italia code **AR22PM**
- Laboratory spray turbine
Speed **300,000 rpm**
FG bur locking system