



## Amigel hand gel disinfectant, 250 ml bottle and dispenser

AMIGEL hand sanitising gel is a ready-to-use, no-rinse solution in a 250ml pump bottle for quick and practical hand hygiene.

**Available**

**SKU:** DIS0250G

**Categories:** Professional Equipment

### Description

---

AMIGEL hand sanitising gel is a ready-to-use solution developed for thorough cleansing of the hands and forearms. Its practical formula is designed for fast, effective application, making it suitable for professional environments where hygiene, ease of use and frequent application are essential.

When used with an automatic wall-mounted dispenser, AMIGEL helps reduce the risk of cross-contamination and supports safer hygiene management in shared and high-traffic settings.

Apply 3ml to dry hands and rub thoroughly for 30 seconds until fully distributed. No rinsing is required after use.

#### Key features AMIGEL hand sanitising gel

- Ready to use
- Fast and practical application
- Suitable for hands and forearms
- No water required
- Ideal for frequent use
- Suitable for use with automatic wall-mounted dispensers
- Fragrance-free
- Helps maintain the skin's natural pH balance
- Formulated with high-grade alcohol to help prevent skin dryness