


Arm, articulated, 400 mm in length,dental white



Monitor support arm for dental unit application, with joint and 35mm spigot, 400mm length.

Available

SKU: BR409010

Categories: Extension Arms

Description

BR409010: articulated monitor support arm with joint and 35mm spigot, 400mm length.

Main features

Vertical pole mounting
Internal routing for electrical cables

Technical specifications

Arm spigot diameter ?35mm
Dimensions L x H: 400 x 190mm
Material: metal
Colour: dental white
Certification: 2001/95/EC
Warranty: 36 months

Application areas

Dental / podiatric treatment units
Equipped poles

Configuration note

The monitor support arm can be combined with VESA mounts, blocks or jaws compatible with a 35mm spigot.