




Printed on 02/06/2025

<https://shop.tecnomeditalia.com/en/product/arm-for-monitor-for-42-45-mm-polecomplete-dental-white/>



Arm for monitor, for 42/45 mm pole,complete, dental white

Complete monitor arm for 42/45 mm pole, color dental white, 330 mm long.

121,19 €

Available

SKU: BR45EC

Categories: Professional Equipment

Description

Monitor arm for 42/45 mm pole

Though they are indispensable items of equipment in any modern dental studio, monitors and computers must not compromise the high levels of hygiene needed. Medical PCs and monitors must therefore conform to strict hygiene standards as well as providing a high level of functionality in terms of mobility and usability.

Stylish, attractive looks are also a plus. Our arms support all types of monitor. Installation is extremely easy and does not require specialist assistance, just a simple Allen key! They also have a convenient handle so that you can move the monitor or all in one PC around easily.

Fixing holes are provided for Vesa 75 and Vesa 100 mountings. This new range of modular monitor arms from Tecnomed Italia offers amazing versatility and lets you specify components of various sizes. You can therefore configure your monitor arm exactly as you want it. Cables can be routed through the arm as an option. Various cables are available, including VGA, USB and power

Technical features Monitor arm for 42/45 mm pole.

Color	dental white
L	330 mm
Yoke ?	42-45 mm
Hollow arm	cable duct
Hollow yoke	no
VESA bracket	75-100
Warranty	24 months