



## Articulated arm for OKKIO camera with base

3 joints articulated arm for OKKIO camera with round base, compatible with standard commercial cameras weighing up to 500 grams.

**Shipping expected in approximately 10 days**

**SKU:** OF933BA

**Categories:** Configurations

### Description

---

**3-joint articulated arm with metal base for the OKKIO medical camera**, designed to ensure **stable, precise and functional positioning** of the camera during dental procedures.

The system allows the **OKKIO medical camera** to be oriented following the clinician's natural line of sight, making it easy to capture **photos and videos** without interfering with clinical workflow.

Thanks to its articulated structure, the camera can be positioned flexibly and maintained at the correct distance from the operating area, ensuring **sharp images and stable recordings**. The base provides a **practical and reliable solution** for clinical documentation, allowing professionals to work **independently**, without interruptions and without the need for a dedicated operator.

An ideal choice for keeping the workstation **clean, professional and well organized**, improving daily efficiency and the overall quality of visual documentation.

The articulated arm can be configured according to operational needs and room characteristics. Thanks to a **dedicated adapter (not included)**, it is also compatible with other cameras featuring a **standard photographic mount**, supporting a maximum weight of **500g**.