


Bearing, high speed, H2,87×3,175×7,520 phen.,steel,groove,MIDWEST 361M



Bearing for high speed handpiece 3.175×7.520×2.870mm, stepped, Tecnomed Italia, with phenolic cage and ceramic balls, compatible with Midwest Stylus 540/541.

Available

SKU: CU36M

Categories: 3,175x7,520x2,87 groove

Description

Stepped bearing for high speed handpiece 3.175×7.520×2.870mm, Tecnomed Italia, spare part for dental turbine, with phenolic cage and ceramic balls.

High speed handpiece compatibility

- Midwest Stylus 540 / 541: bur side and push side

Technical features

- Dimensions: 3.175×7.520×2.870mm
- Stepped bearing
- Phenolic cage
- Ceramic balls
- Radial contact
- ABEC 7 class
- Timken cross reference DA47Z4G-801