


Bearing, low speed, 10x19x7 steel, ZZ



Bearing for low speed handpiece micromotor, 10x19x7mm dimensions, smooth, steel cage, steel balls, Tecnomed Italia brand, bearing closed on both sides, compatible with KaVo K4-K5.

Available

SKU: CM11SH

Categories: Internal diameter 6-12

Description

Low speed handpiece micromotor bearing 10x19x7mm, smooth, by Tecnomed Italia, spare bearing for laboratory micromotors, with steel cage and steel balls.

Laboratory micromotor compatibility

- KaVo K4 – K5

Technical features

- Dimensions: 10x19x7mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Bearing closed on both sides
- Tecnomed Italia brand
- 12-month warranty