



## Low speed handpiece bearing 2.35x5.5x2mm, smooth, steel balls, 200000 rpm

Bearing for low speed handpiece, dimensions 2.35x5.5x2mm, smooth, steel cage, steel balls, closed on one side, 200.000 rpm, Tecnomed Italia brand.

**Available**

**SKU:** CM202/2

**Categories:** Internal diameter 2-3

### Description

---

Low speed handpiece bearing, 2.35x5.5x2mm, smooth type, Tecnomed Italia, spare part for dental handpieces, with steel cage and steel balls.

### Technical features

- Dimensions: height 2mm, internal diameter 2.35mm, external diameter 5.5mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- bearing closed on one side
- 200.000 rpm
- Tecnomed Italia brand
- 12-month warranty