

## Low speed handpiece bearing 2.35x5x1.5mm, smooth, steel balls



Bearing for low speed handpiece, dimensions 2.35x5x1.5mm, smooth, steel cage, steel balls, bearing closed on one side, Tecnomed Italia brand.

**Available**

**SKU:** CM201

**Categories:** Internal diameter 2-3

### Description

---

Low speed handpiece bearing, 2.35x5x1.5mm, smooth type, Tecnomed Italia, spare part for dental handpieces, with steel cage and steel balls.

### Technical features

- Dimensions: height 1.5mm, internal diameter 2.35mm, external diameter 5mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- bearing closed on one side
- Tecnomed Italia brand
- 12-month warranty