

## Low speed handpiece bearing 2.38x4.736x1.588mm, smooth, steel balls



Bearing for low speed handpiece, dimensions 2.38x4.736x1.588mm, smooth, steel cage, steel balls, closed on one side, Tecnomed Italia brand.

**Shipping expected in approximately 15 days**

**SKU:** CM227R

**Categories:** Internal diameter 2-3

### Description

---

Low speed handpiece bearing, 2.38x4.736x1.588mm, smooth type, Tecnomed Italia, spare part for dental handpieces, with steel cage and steel balls.

### Technical features

- Dimensions: height 1.588mm, internal diameter 2.38mm, external diameter 4.736mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Bearing closed on one side
- Tecnomed Italia brand
- 12-month warranty