

Bearing, low speed, 2,380×4,736×2,380 steel



Available

SKU: CM227

Categories: Internal diameter 2-3

Description

Low speed handpiece bearing, 2.38×4.736×1.588mm, Tecnomed Italia smooth type, is a spare part for dental handpieces, with steel cage and steel balls.

Handpiece compatibility

- Ivory SL116C 16:1 Contra-angle Green Ring

Technical features

- Dimensions: height 1.588mm, internal diameter 2.38mm, external diameter 4.736mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- Tecnomed Italia brand
- 12-month warranty