

Low speed handpiece bearing 2.38×4.736×2.38mm, smooth, steel balls



Bearing for low speed handpiece, dimensions 2.38×4.736×2.38mm, smooth, steel cage, steel balls, Tecnomed Italia brand, closed on both sides, compatible with Ivory.

Available

SKU: CM238

Categories: Internal diameter 2-3

Description

Low speed handpiece bearing, 2.38×4.736×2.38mm, Tecnomed Italia smooth type, is a spare part for dental handpieces, with steel cage and steel balls.

Handpiece compatibility

- Ivory SL4201C 20:1 Contra-angle Green Ring

Technical features

- Dimensions: height 2.38mm, internal diameter 2.38mm, external diameter 4.736mm
- Smooth bearing
- Steel cage
- Steel balls
- Radial contact
- Bearing closed on both sides
- Tecnomed Italia brand
- 12-month warranty