


Low speed handpiece bearing 2x6x2.3mm, smooth flanged, steel balls



Bearing for low speed handpiece, dimensions 2x6x2.3mm, smooth flanged, steel cage, steel balls, closed on both sides, Tecnomed Italia brand.

Shipping expected in approximately 15 days

SKU: CM200

Categories: Internal diameter 2-3

Description

Low speed handpiece bearing, 2x6x2.3mm, smooth flanged type, Tecnomed Italia, spare part for dental handpieces, with steel cage and steel balls.

Technical features

- Dimensions: height 2.3mm, internal diameter 2mm, external diameter 6mm
- Smooth flanged bearing
- Steel cage
- Steel balls
- Radial contact
- ABEC 7 class
- Bearing closed on both sides
- Tecnomed Italia brand
- 12-month warranty